

## Phospho-EGFR (Thr678) Rabbit mAb


ABZ4450

### Product Overview

|                             |                                              |
|-----------------------------|----------------------------------------------|
| <b>Name</b>                 | Phospho-EGFR (Thr678) Rabbit mAb             |
| <b>Accession(Primary)</b>   | P00533                                       |
| <b>Host Species</b>         | Rabbit                                       |
| <b>Clonality</b>            | Monoclonal                                   |
| <b>Isotype</b>              | IgG                                          |
| <b>Species Reactivity</b>   | Human                                        |
| <b>Calculated MW</b>        | Calculated MW: 134 kDa; Observed MW: 175 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:1000                              |

### Target information (P00533)

|                          |                                                                                                                                                                                                                                                    |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | EGFR                                                                                                                                                                                                                                               |
| <b>About Gene Symbol</b> | EGFR (Epidermal Growth Factor Receptor) is a receptor tyrosine kinase involved in immune signaling. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Western blot analysis of Phospho-EGFR (Thr678)  
in HeLa lysates using Phospho-EGFR (Thr678)  
antibody.