

Phospho-NFAT2 (Ser237) Rabbit mAb

ABZ4454


Product Overview

Name	Phospho-NFAT2 (Ser237) Rabbit mAb
Accession(Primary)	O95644
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 101 kDa; Observed MW: 101 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:100-1:200

Target information (O95644)

Gene Name NFATC1

About Gene Symbol



Western blot analysis of NFAT2 (Phospho-Ser237) in K562, 3T3 lysates using Phospho-NFAT2 (Ser237) antibody.

