

Phospho-MSK1 (Ser360) Rabbit mAb

ABZ4456


Product Overview

Name	Phospho-MSK1 (Ser360) Rabbit mAb
Accession(Primary)	O75582
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 90 kDa; Observed MW: 90 kDa
Recommended Dilution	WB-1:2000-1:5000; ICC/IF-1:100-1:200


Target information (O75582)

Gene Name RPS6KA5

About Gene Symbol



Western blot analysis of Phospho-MSK1 (Ser360) in C6, CHO-K1 lysates using Phospho-MSK1 (Ser360) antibody.



Immunocytochemistry analysis of MSK1 (Phospho-S360)(green) in Hela using MSK1 (Phospho- S360) antibody, and DAPI(blue)