

Phospho-IRE1 (Ser724) Rabbit mAb

ABZ4459


Product Overview

Name	Phospho-IRE1 (Ser724) Rabbit mAb
Accession(Primary)	O75460
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 110 kDa; Observed MW: 110 kDa
Recommended Dilution	WB-1:1000-1:10000; IP-1:20-1:50

Target information (O75460)

Gene Name ERN1

About Gene Symbol



Western blot analysis of Phospho-IRE1 (Ser724) in MCF-7, HL-60 lysates using Phospho-IRE1 (Ser724) antibody.

