

Phospho-JAK2 (Tyr1007) Rabbit pAb


ABZ4460

Product Overview

Name	Phospho-JAK2 (Tyr1007) Rabbit pAb
Accession(Primary)	O60674
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 131 kDa; Observed MW: 131 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:20000

Target information (O60674)

Gene Name	JAK2
About Gene Symbol	



Western blot analysis of Phospho-JAK2 (Tyr1007) in NIH3T3 lysates using Phospho-JAK2 (Tyr1007) antibody.