

Phospho-FOXO3A (Ser253) Rabbit mAb


ABZ4461

Product Overview


Name	Phospho-FOXO3A (Ser253) Rabbit mAb
Accession(Primary)	O43524
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 71 kDa; Observed MW: 82 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100

Target information (O43524)

Gene Name	FOXO3
About Gene Symbol	FOXO3 (Forkhead Box O3) is a transcription factor involved in gene transcription and regulation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Phospho-FOXO3A (Ser253) in 3T3 lysates using Phospho-FOXO3A (Ser253) antibody.



Western blot analysis of Phospho-FOXO3A (Ser253) in K562, C6 lysates using Phospho-FOXO3A (Ser253) antibody.