

Phospho-TAK1 (Thr187) Rabbit pAb


ABZ4462

Product Overview

Name	Phospho-TAK1 (Thr187) Rabbit pAb
Accession(Primary)	O43318
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 67 kDa; Observed MW: 60 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:10000

Target information (O43318)

Gene Name	MAP3K7
About Gene Symbol	MAP3K7 (Mitogen-Activated Protein Kinase Kinase Kinase 7) is a mitogen-activated protein kinase involved in MAPK/ERK signaling and TGF-β signaling. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Phospho-TAK1 (Thr187)
in 293T lysates using Phospho-TAK1 (Thr187)
antibody.