

## Phospho-TAK1 (Ser439) Rabbit mAb


ABZ4463

### Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | Phospho-TAK1 (Ser439) Rabbit mAb           |
| <b>Accession(Primary)</b>   | O43318                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                          |
| <b>Calculated MW</b>        | Calculated MW: 67 kDa; Observed MW: 78 kDa |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IP-1:10-1:100            |

### Target information (O43318)

|                          |                                                                                                                                                                                                                                                                                                                                                                    |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | MAP3K7                                                                                                                                                                                                                                                                                                                                                             |
| <b>About Gene Symbol</b> | <p>MAP3K7 (Mitogen-Activated Protein Kinase Kinase Kinase 7) is a mitogen-activated protein kinase involved in MAPK/ERK signaling and TGF-<math>\beta</math> signaling. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.</p> |



Western blot analysis of TAK1 (Phospho-Ser439) in 3T3, C6, MCF-7 lysates using Phospho-TAK1 (Ser439) antibody.