

Phospho-IKK alpha/beta (Ser176/Ser177) Rabbit pAb


ABZ4465

Product Overview

Name	Phospho-IKK alpha/beta (Ser176/Ser177) Rabbit pAb
Accession(Primary)	O15111
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 85 kDa; Observed MW: 85 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:10000

Target information (O15111)

Gene Name	CHUK
About Gene Symbol	CHUK (Component Of Inhibitor Of Nuclear Factor Kappa B Kinase Complex) is a protein kinase and signaling regulator involved in NF- κ B signaling. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely st



Western blot analysis of Phospho-IKK alpha/beta (Ser176/Ser177) in 3T3 lysates using Phospho-IKK alpha/beta (Ser176/Ser177) antibody.