

Phospho-c-Jun/JunD (Ser73/100) Rabbit pAb


ABZ4468

Product Overview


Name	Phospho-c-Jun/JunD (Ser73/100) Rabbit pAb
Accession(Primary)	P05412
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 36 kDa; Observed MW: 36 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:10000

Target information (P05412)

Gene Name	JUN/JUND
About Gene Symbol	JUN (Jun Proto-Oncogene, AP-1 Transcription Factor Subunit) is a proto-oncogene involved in gene transcription and regulation. It is an important research target in cancer research and infectious disease and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Phospho-JunD
(Ser73/100) in various lysates using Phospho-JunD
(Ser73/100) antibody.



Western blot analysis of Phospho-JunD
(Ser73/100) in HuvEc lysates using Phospho-
AP1/Jun D (S73/100) antibody.