

Phospho-IGF1 Receptor (Tyr1165/Tyr1166) Rabbit pAb

ABZ4469

Product Overview

Name	Phospho-IGF1 Receptor (Tyr1165/Tyr1166) Rabbit pAb
Accession(Primary)	P08069
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat, Cat
Calculated MW	Calculated MW: 155 kDa; Observed MW: 155 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:20000

Target information (P08069)

Gene Name	IGF1R
About Gene Symbol	

(kD)

