

Phospho-MEK2 (Thr394) Rabbit pAb


ABZ4470

Product Overview


Name	Phospho-MEK2 (Thr394) Rabbit pAb
Accession(Primary)	P36507
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 44 kDa; Observed MW: 44 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; ; ELISA-1:10000; IF-1:50-1:200

Target information (P36507)

Gene Name	MAP2K2
About Gene Symbol	MAP2K2 (Mitogen-Activated Protein Kinase Kinase 2) is a mitogen-activated protein kinase involved in MAPK/ERK signaling. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Phospho-MEK2 (Thr394) in various lysates using Phospho-MEK2 (Thr394) antibody.



Western blot analysis of Phospho-MEK2 (Thr394) in VEC lysates using Phospho-MEK2 (T394) antibody.