

Phospho-JAK2/3 (Tyr966/939) Rabbit pAb


ABZ4471

Product Overview


Name	Phospho-JAK2/3 (Tyr966/939) Rabbit pAb
Accession(Primary)	O60674
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 131 kDa; Observed MW: 120-130 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:50-1:300

Target information (O60674)

Gene Name	JAK2/JAK3
About Gene Symbol	



Western blot analysis of Phospho-JAK2/3 (Tyr966/939) in mouse liver, hela, mouse brain lysates using Phospho-JAK2/3 (Tyr966/939) antibody.



Western blot analysis of Phospho-JAK2/3 (Tyr966/939) in hela lysates using Phospho-JAK2/3 (Tyr966/939) antibody.