

Phospho-GRLF1 (Tyr1087) Rabbit pAb

ABZ4472


Product Overview

Name	Phospho-GRLF1 (Tyr1087) Rabbit pAb
Accession(Primary)	Q9NRY4
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 171 kDa; Observed MW: 190 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:10000

Target information (Q9NRY4)

Gene Name ARHGAP35

About Gene Symbol



Western blot analysis of Phospho-GRLF1 (Tyr1087) in 293T lysates using Phospho-GRLF1 (Tyr1087) antibody.

