

Phospho-Rb (Ser780) Rabbit pAb


ABZ4475

Product Overview


Name	Phospho-Rb (Ser780) Rabbit pAb
Accession(Primary)	P06400
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 106 kDa; Observed MW: 106 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:10000

Target information (P06400)


Gene Name	RB1
About Gene Symbol	RB1 (RB Transcriptional Corepressor 1) is a protein phosphatase involved in gene transcription and regulation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Western blot analysis of Phospho-Rb (Ser780) in various lysates using Phospho-Rb (S780) antibody.



Western blot analysis of Phospho-Rb (Ser780) in KB lysates using Phospho-Rb (S780) antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast carcinoma using Retinoblastoma (Phospho-Ser78) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Sample with blocking peptide on the right.



Western blot analysis of Phospho-Rb (Ser780) in K562 lysates treated with serum 10% using Phospho-Rb (Ser780) antibody. The lane on the right is blocked with the Phospho- peptide.