

Phospho-Sp1 (Thr739) Rabbit pAb


ABZ4477

Product Overview

Name	Phospho-Sp1 (Thr739) Rabbit pAb
Accession(Primary)	P08047
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat, Cat
Calculated MW	Calculated MW: 81 kDa; Observed MW: 81 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:5000

Target information (P08047)

Gene Name	SP1
About Gene Symbol	SP1 (Sp1 Transcription Factor) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Phospho-Sp1 (Thr739) in COS7 lysates using Phospho-Sp1 (Thr739) antibody.