

Phospho-kappa Opioid Receptor (Ser369) Rabbit pAb


ABZ4478

Product Overview


Name	Phospho-kappa Opioid Receptor (Ser369) Rabbit pAb
Accession(Primary)	P41145
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Mouse, Rat
Calculated MW	Calculated MW: 43 kDa; Observed MW: 43 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:10000

Target information (P41145)


Gene Name	OPRK1
About Gene Symbol	




Western blot analysis of Phospho-kappa Opioid Receptor (Ser369) in various lysates using Phospho-kappa Opioid Receptor (Ser369) antibody.



Western blot analysis of Phospho-kappa Opioid Receptor (Ser369) in rat kidney lysates using Phospho-KOR1 (S369) antibody.



Immunohistochemistry analysis of paraffin-embedded Human brain using Phospho-kappa Opioid Receptor (Ser369) antibody. High-pressure and temperature Tris-EDTA pH 8.0 was used for antigen retrieval. Sample with blocking peptide on the right.



Western blot analysis of Phospho-kappa Opioid Receptor (Ser369) in NIH/3T3 lysates using Phospho-kappa Opioid Receptor (Ser369) antibody. The lane on the right is blocked with the Phospho- peptide.