

Phospho-Sodium Potassium ATPase (Ser23) Rabbit pAb


ABZ4479

Product Overview


Name	Phospho-Sodium Potassium ATPase (Ser23) Rabbit pAb
Accession(Primary)	P05023
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Rat
Calculated MW	Calculated MW: 113 kDa; Observed MW: 113 kDa
Recommended Dilution	WB-1:500-1:2000; IF-1:100-1:300

Target information (P05023)


Gene Name	Atp1A1
About Gene Symbol	



Western blot analysis of Phospho-alpha 1 Sodium Potassium ATPase (Ser23) in rat lysates using Phospho-alpha 1 Sodium Potassium ATPase (Ser23) antibody.



Immunofluorescence analysis of Phospho-alpha 1 Sodium Potassium ATPase (Ser23) in NIH/3T3 using ATP1 alpha1/Na+K+ ATPase1 (Phospho-Ser23) antibody. The picture on the right is blocked using the Phospho- peptide.



Western blot analysis of Phospho-alpha 1 Sodium Potassium ATPase (Ser23) in rat brain lysates using ATP1 alpha1/Na+K+ ATPase1 (Phospho-Ser23) antibody. The lane on the right is blocked with the Phospho- peptide.