

# Phospho-Cardiac Troponin I (Ser22/Ser23) Rabbit pAb


ABZ4480

## Product Overview


|                             |                                                                   |
|-----------------------------|-------------------------------------------------------------------|
| <b>Name</b>                 | Phospho-Cardiac Troponin I (Ser22/Ser23) Rabbit pAb               |
| <b>Accession(Primary)</b>   | P19429                                                            |
| <b>Host Species</b>         | Rabbit                                                            |
| <b>Clonality</b>            | Polyclonal                                                        |
| <b>Isotype</b>              | IgG                                                               |
| <b>Species Reactivity</b>   | Mouse, Rat                                                        |
| <b>Calculated MW</b>        | Calculated MW: 24 kDa; Observed MW: 28 kDa                        |
| <b>Recommended Dilution</b> | WB-1:500-1:2000; IHC-1:100-1:300;<br>IF-1:50-1:200; ELISA-1:20000 |

## Target information (P19429)


|                          |       |
|--------------------------|-------|
| <b>Gene Name</b>         | Tnni3 |
| <b>About Gene Symbol</b> |       |




Western blot analysis of Phospho-Cardiac Troponin I (Ser22/Ser23) in SH-SY5Y lysates using Phospho-Cardiac Troponin I (Ser22/Ser23) antibody.



EnzymeLinked Immunosorbent Assay (Phospho-ELISA) for Immunogen Phospho-peptide (Phospho-left) and NonPhospho-peptide (Phospho-right), using TNNI3 (Phospho-Ser22+Ser23) antibody.



Immunohistochemistry analysis of paraffin-embedded Human skeletal muscle using TNNI3 (Phospho-Ser22+Ser23) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Sample with blocking peptide on the right.



Western blot analysis of Phospho-Cardiac Troponin I (Ser22/Ser23) in mouse heart lysates using Phospho-Cardiac Troponin I (Ser22/Ser23) antibody. The lane on the right is blocked with the Phospho-peptide.