

Phospho-p107 (Thr369) Rabbit pAb

ABZ4482


Product Overview

Name	Phospho-p107 (Thr369) Rabbit pAb
Accession(Primary)	P28749
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:5000


Target information (P28749)

Gene Name RBL1

About Gene Symbol



EnzymeLinked Immunosorbent Assay (Phospho-ELISA) for Immunogen Phospho-peptide (Phospho-left) and NonPhospho-peptide (Phospho-right), using RBL1 (Phospho-Thr36antibody



Immunohistochemical analysis of paraffin-embedded Human tonsils using Phospho-p107 (Thr369) antibody. Sample with blocking peptide on the right. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.