

## Phospho-CD31 (Tyr713) Rabbit pAb


ABZ4483

### Product Overview


|                             |                                                                     |
|-----------------------------|---------------------------------------------------------------------|
| <b>Name</b>                 | Phospho-CD31 (Tyr713) Rabbit pAb                                    |
| <b>Accession(Primary)</b>   | P16284                                                              |
| <b>Host Species</b>         | Rabbit                                                              |
| <b>Clonality</b>            | Polyclonal                                                          |
| <b>Isotype</b>              | IgG                                                                 |
| <b>Species Reactivity</b>   | Human, Mouse                                                        |
| <b>Calculated MW</b>        | Calculated MW: 83 kDa; Observed MW: 150 kDa                         |
| <b>Recommended Dilution</b> | WB-1:500-1:2000; IHC-1:100-1:300;<br>IF-1:200-1:1000; ELISA-1:20000 |

### Target information (P16284)

|                          |        |
|--------------------------|--------|
| <b>Gene Name</b>         | PECAM1 |
| <b>About Gene Symbol</b> |        |




EnzymeLinked Immunosorbent Assay (Phospho-ELISA) for Immunogen Phospho-peptide (Phospho-left) and NonPhospho-peptide (Phospho-right), using PECAM1 (Phospho-Tyr71)antibody



HepG2 HepG2

Immunofluorescence analysis of Phospho-CD31 (Tyr713) in HeLa using Phospho-CD31 (Tyr713) antibody. The picture on the right is blocked with the Phospho- peptide.



Western blot analysis of Phospho-CD31 (Tyr713) in HepG2 lysates using Phospho-CD31 (Tyr713) antibody. The lane on the right is blocked with the Phospho- peptide.