

Phospho-Histone H4 (Ser1) Rabbit pAb

ABZ4486

Product Overview

Name	Phospho-Histone H4 (Ser1) Rabbit pAb
Accession(Primary)	P62805
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 11 kDa; Observed MW: 17 kDa
Recommended Dilution	WB-1:1000-1:2000

Target information (P62805)

Gene Name H4C1

About Gene Symbol



