

Phospho-NMDAR2A/2B (Tyr1246/Tyr1252) Rabbit pAb


ABZ4487

Product Overview

Name	Phospho-NMDAR2A/2B (Tyr1246/Tyr1252) Rabbit pAb
Accession(Primary)	Q12879
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:10000

Target information (Q12879)

Gene Name	GRIN2A/GRIN2B
About Gene Symbol	



Immunohistochemistry analysis of paraffin-embedded Human brain using Phospho-NMDAR2A/2B (Tyr1246/Tyr1252) antibody. High-pressure and temperature Tris-EDTA pH 8.0 was used for antigen retrieval. Sample with blocking peptide on the right.