

Phospho-CAD (Thr456) Rabbit pAb

ABZ4488


Product Overview

Name	Phospho-CAD (Thr456) Rabbit pAb
Accession(Primary)	P27708
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Recommended Dilution	IHC-1:100-1:300; ELISA-1:5000; IF-1:50-1:200

Target information (P27708)

Gene Name CAD

About Gene Symbol



Immunohistochemistry analysis of paraffin-embedded Human colon carcinoma tissue using Phospho-CAD (Thr456) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Sample with blocking peptide on the right.