

Lamin A/C Rabbit mAb


ABZ4489

Product Overview


Name	Lamin A/C Rabbit mAb
Accession(Primary)	P02545
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 74 kDa; Observed MW: 74,63 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200; FC-1:20; IP-1:10-1:100

Target information (P02545)


Gene Name	LMNA
About Gene Symbol	




Western blot analysis of Lamin A/C in rat heart, rat skeletal muscle lysates using Lamin A/C antibody.



Western blot analysis of Lamin A/C in 3T3, HeLa lysates using Lamin A/C antibody.



Immunocytochemistry analysis of Lamin A/C(green) in HeLa using Lamin A/C antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using Lamin A/C antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.