

Phospho-HDAC5 (Ser498) Rabbit pAb

ABZ4490

Product Overview

Name	Phospho-HDAC5 (Ser498) Rabbit pAb
Accession(Primary)	Q9UQL6
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 122 kDa; Observed MW: 122 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:20000

Target information (Q9UQL6)

Gene Name	HDAC5
About Gene Symbol	

