

Phospho-PKR (Thr446) Rabbit pAb


ABZ4492

Product Overview

Name	Phospho-PKR (Thr446) Rabbit pAb
Accession(Primary)	P19525
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 62 kDa; Observed MW: 68 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:10000

Target information (P19525)

Gene Name	EIF2AK2
About Gene Symbol	EIF2AK2 (eIF-2A protein kinase 2) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Phospho-PKR (Thr446) in K562 lysates using Phospho-PKR (Thr446) antibody.