

Phospho-OXSR1 (Thr185) Rabbit pAb


ABZ4494

Product Overview


Name	Phospho-OXSR1 (Thr185) Rabbit pAb
Accession(Primary)	O95747
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 58 kDa; Observed MW: 65 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:10000

Target information (O95747)

Gene Name	OXSR1
About Gene Symbol	



Western blot analysis of Phospho-OXSR1 (Thr185) in HepG2 lysates using Phospho-OXSR1 (Thr185) antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast carcinoma tissue using OXSR1 (Phospho-Thr185) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Sample with blocking peptide on the right.