

Phospho-Connexin 43 / GJA1 (Tyr265) Rabbit pAb

ABZ4495

Product Overview

Name	Phospho-Connexin 43 / GJA1 (Tyr265) Rabbit pAb
Accession(Primary)	P17302
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 43 kDa; Observed MW: 43 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:10000

Target information (P17302)

Gene Name GJA1

About Gene Symbol



