

Phospho-ATF4 (Ser219) Rabbit pAb

ABZ4496

Product Overview

Name	Phospho-ATF4 (Ser219) Rabbit pAb
Accession(Primary)	P18848
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 39 kDa; Observed MW: 39 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:10000

Target information (P18848)

Gene Name ATF4

About Gene Symbol



