

Phospho-EPHA4 (Tyr596) Rabbit pAb

ABZ4497


Product Overview

Name	Phospho-EPHA4 (Tyr596) Rabbit pAb
Accession(Primary)	P54764
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 110 kDa; Observed MW: 110 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:10000

Target information (P54764)

Gene Name EPHA4

About Gene Symbol



Western blot analysis of Phospho-EPHA4 (Tyr596) in Jurkat lysates using Phospho-EPHA4 (Tyr596) antibody.

