

Phospho-FGFR4 (Tyr642) Rabbit pAb

ABZ4498

Product Overview

Name	Phospho-FGFR4 (Tyr642) Rabbit pAb
Accession(Primary)	P22455
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 88 kDa; Observed MW: 88 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:10000

Target information (P22455)

Gene Name FGFR4**About Gene Symbol**

