

Phospho-Annexin 2 (Ser26) Rabbit pAb

ABZ4499

Product Overview

Name	Phospho-Annexin 2 (Ser26) Rabbit pAb
Accession(Primary)	P07355
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 39 kDa; Observed MW: 39 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:10000

Target information (P07355)

Gene Name ANXA2

About Gene Symbol



