

Phospho-DRP1 (Ser637) Rabbit pAb

ABZ4500

Product Overview

Name	Phospho-DRP1 (Ser637) Rabbit pAb
Accession(Primary)	O00429
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 82 kDa; Observed MW: 82 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:100-1:300; ELISA-1:20000

Target information (O00429)

Gene Name	DNM1L
About Gene Symbol	DNM1L (Dynamin 1 Like) is a dynamin GTPase involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

