

Phospho-DDIT3 (Ser30) Rabbit pAb


ABZ4502

Product Overview

Name	Phospho-DDIT3 (Ser30) Rabbit pAb
Accession(Primary)	P35638
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 19 kDa; Observed MW: 19 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:20000

Target information (P35638)

Gene Name	DDIT3
About Gene Symbol	DDIT3 (DNA Damage Inducible Transcript 3) is an epigenetic and nuclear regulatory protein involved in DNA damage response and repair. It is an important research target in epigenetic regulation, chromatin biology, and cancer research and is widely studied



Western blot analysis of Phospho-CHOP (Ser30) in Jurkat lysates using Phospho-CHOP (Ser30) antibody.