

Phospho-CSK (Ser364) Rabbit pAb


ABZ4503

Product Overview

Name	Phospho-CSK (Ser364) Rabbit pAb
Accession(Primary)	P41240
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 51 kDa; Observed MW: 55 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:10000

Target information (P41240)

Gene Name	CSK
About Gene Symbol	CSK (C-Terminal Src Kinase) is a receptor tyrosine kinase involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Phospho-CSK (Ser364) in HUVEC lysates using Phospho-CSK (Ser364) antibody.