

Phospho-p38 MAPK (Tyr323) Rabbit pAb


ABZ4504

Product Overview


Name	Phospho-p38 MAPK (Tyr323) Rabbit pAb
Accession(Primary)	Q16539
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 41 kDa; Observed MW: 35 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:5000

Target information (Q16539)

Gene Name	MAPK14
About Gene Symbol	MAPK14 (Mitogen-Activated Protein Kinase 14) is a mitogen-activated protein kinase involved in MAPK/ERK signaling and JAK-STAT/cytokine signaling. It is an important research target in inflammatory disease and is widely studied to elucidate cellular mecha



Western blot analysis of Phospho-p38 (Tyr323) in HepG2 lysates using Phospho-p38 (Tyr323) antibody.



Immunohistochemistry analysis of paraffin-embedded rat kidney tissue using Phospho-p38 (Tyr323) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.