


Phospho-Leptin Receptor (Tyr1141) Rabbit pAb

ABZ4505

Product Overview

Name	Phospho-Leptin Receptor (Tyr1141) Rabbit pAb
Accession(Primary)	P48357
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 132 kDa; Observed MW: 132 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:5000

Target information (P48357)

Gene Name LEPR**About Gene Symbol**

Western blot analysis of Phospho-Leptin Receptor (Tyr1141) in K562 lysates using Phospho-Leptin Receptor (Tyr1141) antibody.

