

Phospho-PTEN (Ser380) Rabbit pAb

ABZ4506


Product Overview

Name	Phospho-PTEN (Ser380) Rabbit pAb
Accession(Primary)	P60484
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Recommended Dilution	IHC-1:100-1:300; ELISA-1:10000; IF-1:50-1:200


Target information (P60484)

Gene Name PTEN

About Gene Symbol



Immunohistochemistry analysis of paraffin-embedded Human breast carcinoma tissue using Phospho-PTEN (Ser380) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Sample with blocking peptide on the right.



Immunofluorescence analysis of Phospho-PTEN (Ser380) in HeLa cells using PTEN (Phospho-Ser380) antibody, and DAPI (blue).