


Phospho-S1P1/EDG1 (Thr236) Rabbit pAb

ABZ4507


Product Overview

Name	Phospho-S1P1/EDG1 (Thr236) Rabbit pAb
Accession(Primary)	P21453
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 43 kDa; Observed MW: 42 kDa
Recommended Dilution	WB-1:500-1:2000; IF-1:200-1:1000

Target information (P21453)

Gene Name S1PR1**About Gene Symbol**

Western blot analysis of Phospho-EDG1 (Thr236) in HepG2, mouse brain lysates using Phospho-EDG1 (Thr236) antibody.



Immunofluorescence analysis of Phospho-EDG1 (Thr236) in COS7 using Phospho-EDG1 (Thr236) antibody. The picture on the right is blocked using the Phospho- peptide.