

Phospho-Synapsin I (Ser9) Rabbit mAb

ABZ4510


Product Overview


Name	Phospho-Synapsin I (Ser9) Rabbit mAb
Accession(Primary)	P17600
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 74 kDa; Observed MW: 77 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-1:50-1:100

Target information (P17600)

Gene Name SYN1

About Gene Symbol





Immunohistochemistry analysis of paraffin-embedded Human Brain using Synapsin I (Phospho- S9) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.