

## Phospho-RSK4 (Ser232) Rabbit mAb


ABZ4511

### Product Overview



|                             |                                                          |
|-----------------------------|----------------------------------------------------------|
| <b>Name</b>                 | Phospho-RSK4 (Ser232) Rabbit mAb                         |
| <b>Accession(Primary)</b>   | Q9UK32                                                   |
| <b>Host Species</b>         | Rabbit                                                   |
| <b>Clonality</b>            | Monoclonal                                               |
| <b>Isotype</b>              | IgG                                                      |
| <b>Species Reactivity</b>   | Human, Mouse                                             |
| <b>Calculated MW</b>        | Calculated MW: 84 kDa; Observed MW: 84 kDa               |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:50-1:100;<br>ICC/IF-1:50-1:100 |

### Target information (Q9UK32)

|                          |         |
|--------------------------|---------|
| <b>Gene Name</b>         | RPS6KA6 |
| <b>About Gene Symbol</b> |         |



Immunocytochemistry analysis of Phospho-RSK4 (Ser232) (green) in LNCap using Phospho-RSK4 (Ser232) antibody, and DAPI (blue).



Western blot analysis of Phospho-RSK4 (Ser232) in HeLa lysates using Phospho-RSK4 (Ser232) antibody.