

Phospho-PTP1B (Ser352) Rabbit mAb

ABZ4514


Product Overview

Name	Phospho-PTP1B (Ser352) Rabbit mAb
Accession(Primary)	P18031
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 50 kDa; Observed MW: 50 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:50-1:100


Target information (P18031)

Gene Name PTPN1

About Gene Symbol



Western blot analysis of Phospho-PTP1B (Ser352) in Jurkat lysates using Phospho-PTP1B (Ser352) antibody.



Immunocytochemistry analysis of Phospho-PTP1B
(Ser352)(green) in HeLa using Phospho-PTP1B
(Ser352) antibody, and DAPI(blue)