

## Phospho-Presenilin 2/AD5 (Ser330) Rabbit mAb

ABZ4515


### Product Overview

|                             |                                              |
|-----------------------------|----------------------------------------------|
| <b>Name</b>                 | Phospho-Presenilin 2/AD5 (Ser330) Rabbit mAb |
| <b>Accession(Primary)</b>   | P49810                                       |
| <b>Host Species</b>         | Rabbit                                       |
| <b>Clonality</b>            | Monoclonal                                   |
| <b>Isotype</b>              | IgG                                          |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                            |
| <b>Calculated MW</b>        | Calculated MW: 50 kDa; Observed MW: 23 kDa   |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000                             |

### Target information (P49810)

**Gene Name** PSEN2

**About Gene Symbol**



Western blot analysis of Phospho-Presenilin 2 (Ser330) in K562, C6, 3T3, HeLa lysates using Phospho-Presenilin 2 (Ser330) antibody.

