

Phospho-BRAF (Thr401) Rabbit mAb

ABZ4517

Product Overview

Name	Phospho-BRAF (Thr401) Rabbit mAb
Accession(Primary)	P15056
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 84 kDa; Observed MW: 84 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (P15056)

Gene Name BRAF

About Gene Symbol



