

Phospho-CBL (Tyr774) Rabbit mAb

ABZ4518


Product Overview

Name	Phospho-CBL (Tyr774) Rabbit mAb
Accession(Primary)	P22681
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 100 kDa; Observed MW: 120 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:50-1:100; FC-1:20

Target information (P22681)

Gene Name CBL

About Gene Symbol



Western blot analysis of Phospho-CBL (Tyr774) in K562, C6, HeLa lysates using Phospho-CBL (Tyr774) antibody.

