

Phospho-EIF2S1 (Ser51) Rabbit mAb


ABZ4519

Product Overview

Name	Phospho-EIF2S1 (Ser51) Rabbit mAb
Accession(Primary)	P05198
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 36 kDa; Observed MW: 38-45 kDa
Recommended Dilution	WB-1:1000-1:10000; IHC-1:200-1:800

Target information (P05198)

Gene Name	EIF2S1
About Gene Symbol	EIF2S1 (Eukaryotic Translation Initiation Factor 2 Subunit Alpha) is a protein kinase and signaling regulator involved in translational regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Phospho-eIF2A (Ser51) in K562, C6, 3T3, HeLa lysates using Phospho-eIF2A (Ser51) antibody.