

Phospho-c-Jun (Ser63) Rabbit mAb


ABZ4520

Product Overview


Name	Phospho-c-Jun (Ser63) Rabbit mAb
Accession(Primary)	P05412
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 36 kDa; Observed MW: 43 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200

Target information (P05412)

Gene Name	JUN
About Gene Symbol	JUN (Jun Proto-Oncogene, AP-1 Transcription Factor Subunit) is a proto-oncogene involved in gene transcription and regulation. It is an important research target in cancer research and infectious disease and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Phospho-c-Jun (Ser63) in K562, C6, 3T3 lysates using Phospho-c-Jun (Ser63) antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using c-Jun (Phospho- S63) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.