

Phospho-Tau (Ser202/Thr205) Rabbit mAb

ABZ4522


Product Overview

Name	Phospho-Tau (Ser202/Thr205) Rabbit mAb
Accession(Primary)	P10636
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 79 kDa; Observed MW: 50-80 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20-1:50

Target information (P10636)

Gene Name MAPT

About Gene Symbol



Western blot analysis of Tau (Phospho-Ser202/Thr205) in Human Fetal brain lysates using Phospho-Tau (Ser202/Thr205) antibody.

