


Phospho-Tau (Thr205) Rabbit mAb

ABZ4523

Product Overview

Name	Phospho-Tau (Thr205) Rabbit mAb
Accession(Primary)	P10636
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 79 kDa; Observed MW: 50-80 kDa
Recommended Dilution	WB-1:2000-1:5000

Target information (P10636)

Gene Name MAPT**About Gene Symbol****Human Fetal brain**

Western blot analysis of Tau (Phospho-Thr205) in Human Fetal brain lysates using Phospho-Tau (Thr205) antibody.

